

# EpiQuik™ Global Di-Methyl Histone H4R3 Quantification Kit (Colorimetric)

Base Catalog # P-3090

## PLEASE READ THIS ENTIRE USER GUIDE BEFORE USE

**Uses:** The EpiQuik™ Global Di-Methyl Histone H4R3 Quantification Kit (Colorimetric) is suitable for specifically measuring global histone H4 arginine 3 di-methylation from a broad range of species such as mammals, plants, fungi, and bacteria, in a variety of forms including cultured cells and fresh tissues. Histone extracts can be prepared by using your own successful method. For your convenience and the best results, Epigentek offers a histone extraction kit (Cat. # OP-0006) optimized for use with this kit. Histone extracts can be used immediately or stored at  $-80^{\circ}\text{C}$  for future use.

**Input Material:** Input materials can be histone extracts or nuclear extracts. The amount of histone extracts for each assay can be 0.1  $\mu\text{g}$  to 2  $\mu\text{g}$  with an optimal range of 0.5 to 1  $\mu\text{g}$ .

**Internal Control:** The assay control (methylated histone H4-Arg 3) is provided in this kit for the quantification of global di-methyl histone H4R3. Because content of di-methyl histone H4R3 can vary from tissue to tissue, and from normal and diseased states, it is advised to run replicate samples to ensure that the signal generated is validated.

**Precautions:** To avoid cross-contamination, carefully pipette the sample or solution into the strip wells. Use aerosol-barrier pipette tips and always change pipette tips between liquid transfers. Wear gloves throughout the entire procedure. In case of contact between gloves and sample, change gloves immediately.

## KIT CONTENTS

Component	48 Assays Cat. #P-3090-48	96 Assays Cat. #P-3090-96	Storage Upon Receipt
<b>WB</b> (10X Wash Buffer)	14 ml	28 ml	4°C
<b>HB</b> (Histone Buffer)	4 ml	8 ml	4°C
<b>BB</b> (Blocking Buffer)	10 ml	20 ml	4°C
<b>MER3</b> (Capture Antibody, 1000X)*	5 µl	10 µl	4°C
<b>RDA</b> (Detection Antibody, 2000X)*	6 µl	12 µl	-20°C
<b>ES</b> (Enhancer solution)	6 µl	12 µl	-20°C
<b>DS</b> (Developer Solution)	5 ml	10 ml	4°C
<b>SS</b> (Stop Solution)	5 ml	10 ml	RT
Di-Methyl H4R3 control ( 50 µg/ml)	10 µl	20 µl	-20°C
8-Well Assay Strips (With Frame)	6	12	4°C
Adhesive Covering Film	1	1	RT
User Guide	1	1	RT

\* Spin the solution down to the bottom prior to use.

## SHIPPING & STORAGE

The kit is shipped in three parts: the first part at ambient room temperature, and the second and third parts on frozen ice packs at 4°C. Upon receipt: (1) Store **RDA**, **ES**, and **Di-Methyl H4R3 control** at -20°C away from light; (2) Store **WB**, **HB**, **BB**, **MER3**, **DS**, and **8-Well Assay Strips** at 4°C away from light; and (3) Store remaining components (**SS** and **Adhesive Covering Film**) at room temperature away from light.

All components of the kit are stable for 6 months from the date of shipment, when stored properly.

*Note: (1) Check if **WB** (10X Wash Buffer) contains salt precipitates before use. If so, warm (at room temperature or 37°C) and shake the buffer until the salts are re-dissolved; and (2) check if a blue color present in **DS** (Developer Solution), which would indicate contamination of the solution and should not be used. To avoid contamination, transfer the amount of **DS** required into a secondary container (tube or vial) before adding **DS** into the assay wells.*

## MATERIALS REQUIRED BUT NOT SUPPLIED

- ☐ Adjustable pipette or multiple-channel pipette
- ☐ Multiple-channel pipette reservoirs
- ☐ Aerosol resistant pipette tips
- ☐ Microplate reader capable of reading absorbance at 450 nm
- ☐ 1.5 ml microcentrifuge tubes
- ☐ Incubator for 37°C incubation

- ☐ Distilled water
- ☐ Histone extracts
- ☐ Parafilm M or aluminum foil

## GENERAL PRODUCT INFORMATION

**Quality Control:** Each lot of the EpiQuik™ Global Di-Methyl Histone H4R3 Quantification Kit (Colorimetric) is tested against predetermined specifications to ensure consistent product quality. Epigentek guarantees the performance of all products in the manner described in our product instructions.

**Product Warranty:** If this product does not meet your expectations, simply contact our technical support unit or your regional distributor. We also encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

**Safety:** Suitable lab coat, disposable gloves, and proper eye protection are required when working with this product.

**Product Updates:** Epigentek reserves the right to change or modify any product to enhance its performance and design. The information in this User Guide is subject to change at any time without notice. Thus, only use the User Guide that was supplied with the kit when using that kit.

**Usage Limitation:** The EpiQuik™ Global Di-Methyl Histone H4R3 Quantification Kit (Colorimetric) is for research use only and is not intended for diagnostic or therapeutic application.

## A BRIEF OVERVIEW

Arginine histone methylation is one of the many important epigenetic marks, and is essential for the regulation of multiple cellular processes. Arginine methylation of histones H3 (Arg2, 8, 17, 26) and H4 (Arg3) promotes transcriptional activation and is mediated by a family of protein arginine methyltransferases (PRMTs). There are 9 types of PRMTs found in humans but only 7 members are reported to methylate histones. They can mediate mono or dimethylation of arginine residues. These enzymes use S-adenosyl-methionine (SAM) as a methyl donor and transfer it to the guanidinium side chain of arginine. Based on the position of methyl group addition, the PRMTs can be classified into type I (CARM1, PRMT1, PRMT2, PRMT3, PRMT6, and PRMT8) and type II (PRMT5 and PRMT7).






Fig. 1. Histone arginine methylation reaction catalyzed by PRMTs.


Symmetric di-methylation of histone H4 arg3 (H4R3) is catalyzed by type II PRMTs, which are found to be strongly implicated in diseases like cancer. For example, PRMT5 plays a role in the repression of certain tumor suppressor genes such as RB tumor suppressors while PRMT7 overexpression is observed in breast cancer. The global H4R3 di-methylation can be changed by inhibition or activation of type II PRMTs. Therefore, quantitative detection of global symmetric di-methyl histone H4R3 would provide useful information for better understanding epigenetic regulation of gene activation and silencing, as well as for developing PRMT-targeted drugs.

The EpiQuik™ Global Di-Methyl Histone H4R3 Quantification Kit (Colorimetric) is designed to quantitatively detect global di-methyl histone H4R3. This kit has the following advantages:


- Quick and efficient procedure, which can be finished within 3.5 hours.
- Innovative colorimetric assay without the need for radioactivity, electrophoresis, or chromatography.
- Specifically captures symmetric di-methylated H4R3 with the detection limit as low as 0.5 ng/well and detection range from 100 ng to 2 µg/well of histone extracts.
- The control is conveniently included for the quantification of di-methylated H4R3.
- Strip microplate format makes the assay flexible: manual or high throughput.
- Simple, reliable, and consistent assay conditions.

## PRINCIPLE & PROCEDURE


The EpiQuik™ Global Di-Methyl Histone H4R3 Quantification Kit (Colorimetric) is designed for measuring global histone H4R3 di-methylation. In an assay with this kit, the histone proteins are stably spotted on the strip wells. The di-methyl histone H4R3 can be recognized with a high-affinity antibody and detected with a detection antibody, followed by a color development reagent. The ratio of di-methylated H4R3 is proportional to the intensity of absorbance. The absolute amount of di-methylated H4R3 can be quantitated by comparing to the standard control.



Schematic procedure of the EpiQuik™ Global Di-Methyl Histone H4R3 Quantification Kit (Colorimetric)



Histone extracts were prepared from MDA-231 cells using the EpiQuik™ Total Histone Extraction Kit and the amount of dimethyl-H4R3 was measured using the EpiQuik™ Global Dimethyl Histone H4R3 Quantification Kit (Colorimetric).



Illustrated standard curve generated with di-methyl H4R3 control.

## PROTOCOL

For the best results, please read the protocol in its entirety prior to starting your experiment

### Starting Materials

**Input Amount:** The amount of histone extracts for each assay can be between 0.1 µg and 2 µg with an optimal range of 0.2 to 0.5 µg.



*Histone Extraction:* You can use your method of choice for preparing histone extracts from the treated and untreated samples. Epigentek also offers a histone extraction kit (Cat # OP-0006) optimized for use with this kit.

Histone extracts should be stored in aliquots at  $-80^{\circ}\text{C}$  until use.

### **1. Working Buffer and Solution Preparation**

- a. Prepare **Diluted WB** 1X Wash Buffer:

48-Assay Kit: Add 13 ml of **WB** 10X Wash Buffer to 117 ml of distilled water and adjust pH to 7.2-7.5.

96-Assay Kit: Add 26 ml of **WB** 10X Wash Buffer to 234 ml of distilled water and adjust pH to 7.2-7.5.

This **Diluted WB** 1X Wash Buffer can now be stored at  $4^{\circ}\text{C}$  for up to six months.

- b. Prepare **Diluted MER3** Capture Antibody Solution:

Dilute **MER3** Capture Antibody with **Diluted WB** 1X Wash Buffer at a ratio of 1:1000 (i.e., add 1  $\mu\text{l}$  of **MER3** to 1000  $\mu\text{l}$  of **Diluted WB**). 50  $\mu\text{l}$  of **Diluted MER3** will be required for each assay well.

- c. Prepare **Diluted RDA** Detection Antibody Solution:

Dilute **RDA** Detection Antibody with **Diluted WB** 1X Wash Buffer at a ratio of 1:2000 (i.e., add 1  $\mu\text{l}$  of **RDA** to 2000  $\mu\text{l}$  of **Diluted WB**). 50  $\mu\text{l}$  of **Diluted RDA** will be required for each assay well.

- d. Prepare **Diluted ES** Enhancer Solution:

Dilute **ES** Enhancer Solution with **Diluted WB** 1X Wash Buffer at a ratio of 1:5000 (i.e., add 1  $\mu\text{l}$  of Wash Buffer at a ratio of 1:5000 (i.e., add 1  $\mu\text{l}$  of **ES** to 5000  $\mu\text{l}$  of **WB**). About 50  $\mu\text{l}$  of this **Diluted ES** will be required for each assay well.

- e. Prepare **Diluted di-methyl H4R3 Control** Standard

Suggested Standard Curve Preparation: First, dilute **H4R3 Control** with **HB** histone buffer to 10 ng/ $\mu\text{l}$  by adding 2  $\mu\text{l}$  of **H4R3 Control** to 8  $\mu\text{l}$  of **HB** histone buffer. Then, further prepare five concentrations by combining the 10 ng/ $\mu\text{l}$  **Diluted H4R3 Control** with **HB** into final concentrations of 0.5, 1, 2, 5, and 10 ng/ $\mu\text{l}$  according to the following dilution chart:

Tube	H4R3 (20 ng/ $\mu\text{l}$ )	HB	Resulting H4R3 Concentration
1	1.0 $\mu\text{l}$	19.0 $\mu\text{l}$	0.5 ng/ $\mu\text{l}$
2	1.0 $\mu\text{l}$	9.0 $\mu\text{l}$	1 ng/ $\mu\text{l}$
3	1.0 $\mu\text{l}$	4.0 $\mu\text{l}$	2 ng/ $\mu\text{l}$
4	2.0 $\mu\text{l}$	2.0 $\mu\text{l}$	5 ng/ $\mu\text{l}$
5	4.0 $\mu\text{l}$	0.0 $\mu\text{l}$	10 ng/ $\mu\text{l}$

*Note: Keep each of the diluted solutions except **WB** 1X Wash Buffer on ice until use. Any remaining diluted solutions other than **Diluted WB** should be discarded if not used within the same day.*

## **2. Histone Binding**

- a. Predetermine the number of strip wells required for your experiment. It is advised to run replicate samples (include blank and positive controls) to ensure that the signal generated is validated. Carefully remove un-needed strip wells from the plate frame and place them back in the bag (seal the bag tightly and store at 4°C).
- b. **Blank Wells:** Add 49 µl of **HB** to each blank well.
- c. **Standard Wells:** Add 49 µl of **HB** and 1 µl of **Diluted H4R3 control** to each standard well with a minimum of five wells, each at a different concentration between 0.5 and 10 ng/µl (based on the dilution chart in Step 1e; see **Table 2** under the “Suggested Strip Well Setup” section as an example).
- d. **Sample Wells:** Add 46 to 49 µl of **HB** and 1 to 4 µl of your histone extracts. Total volume should be 50 µl per well.

*Note: (1) Follow the suggested well setup diagrams; (2) It is recommended to use 0.2 µg to 0.5 µg of histone extract per well.*

- e. Tightly cover strip-well microplate with **Adhesive Covering Film** to avoid evaporation and incubate at 37°C for 90 to 120 min.

*Note: The **Adhesive Covering Film** can be cut to the required size to cover the strips based on the number of strips to be used.*

- f. Remove the reaction solution from each well. Add 150 µl of **BB** block buffer to each well, then cover with Parafilm M or aluminium foil and incubate at 37°C for 30 min.
- g. Remove the reaction solution from each well. Wash each well three times with 150 µl of the **Diluted WB 1X Wash Buffer** each time.

## **3. Antibody Binding and Signal Enhancing**

- a. Add 50 µl of the **Diluted MER3** to each well, then cover with Parafilm M or aluminium foil and incubate at room temperature for 60 min.
- b. Remove the **Diluted MER3** solution from each well.
- c. Wash each well three times with 150 µl of the **Diluted WB** each time.
- d. Add 50 µl of the **Diluted RDA** to each well, then cover with Parafilm M or aluminium foil and incubate at room temperature for 30 min.
- e. Remove the **Diluted RDA** solution from each well.
- f. Wash each well four times with 150 µl of the **Diluted WB** each time.
- g. Add 50 µl of the **Diluted ES** to each well, then carefully cover with Parafilm M or aluminium foil and incubate at room temperature for 30 min.
- h. Remove the **Diluted ES** solution from each well.
- i. Wash each well with 150 µl of the **Diluted WB** each time for five times.

*Note: Ensure any residual wash buffer in the wells is removed as much as possible at each wash step.*

#### **4. Signal Detection**

- Add 100 µl of **DS** to each well and incubate at room temperature for 1 to 10 min away from light. Begin monitoring color change in the sample wells and control wells. The **DS** solution will turn blue in the presence of sufficient methylated products.
- Add 100 µl of **SS** to each well to stop enzyme reaction when color in the positive control wells turns medium blue. The color will change to yellow after adding **SS** and the absorbance should be read on a microplate reader within 2 to 10 min at 450 nm with an optional reference wavelength of 655 nm.

*Note: (1) Most microplate readers have the capability to carry out dual wavelength analysis and will automatically subtract reference wavelength absorbance from the test wavelength absorbance. If your plate reader does not have this capability, the plate can be read twice, once at 450 nm and once at 655 nm. Then, manually subtract the 655 nm ODs from 450 nm ODs; (2) If the strip-well microplate frame does not fit in the microplate reader, transfer the solution to a standard 96-well microplate.*

#### **5. Di-Methyl H4R3 Calculation**

- Calculate the average duplicate readings for the sample wells and blank wells.
- Calculate % histone H4R3 di-methylation change using the following formula:

$$\text{Di-Methyl H4R3\%} = \frac{\text{Treated (Tested) Sample OD} - \text{Blank OD}}{\text{Untreated (Control) Sample OD} - \text{Blank OD}} \times 100\%$$

Example calculation:

Average OD450 of treated sample is 0.5  
 Average OD450 of untreated control is 0.9  
 Average OD450 of blank is 0.1


$$\text{Di-Methyl H4R3 \%} = \frac{(0.5 - 0.1)}{0.9 - 0.1} \times 100\% = 50\%$$

For accurate calculation:

- Generate a standard curve and plot OD value versus amount of **H4R3 control** at each concentration point.
- Determine the slope as OD/ng (you can use Microsoft Excel statistical functions for slope calculation), then calculate the amount of di-methyl H4R3 using the following formulas:

$$\text{Di-Methyl H4R3 (ng/mg protein)} = \frac{(\text{Sample OD} - \text{Blank OD})}{\text{Slope} \times \text{Protein Amount (ug*)}} \times 1000$$

\* Histone extract added into sample wells at step 2d.





## SUGGESTED BUFFER AND SOLUTION SETUP

**Table 1.** Approximate amount of required buffers and solutions for defined assay wells based on the protocol.

Reagents	1 well	1 strip (8 wells)	2 strips (16 wells)	6 strips (48 wells)	12 strips (96 wells)
Diluted WB	2.5 ml	20 ml	40 ml	120 ml	240 ml
HB	50 µl	400 µl	800 µl	2400 µl	4800 µl
BB	0.15 ml	1.2 ml	2.5 ml	7.5 ml	14.5 ml
H4R3 control	N/A	N/A	4 µL (optional)	8 µl	8 µl
Diluted MER3	50 µl	400 µl	800 µl	2400 µl	4800 µl
Diluted RDA	50 µl	400 µl	800 µl	2400 µl	4800 µl
Diluted ES	50 µl	400 µl	800 µl	2400 µl	4800 µl
Developer Solution	0.1 ml	0.8 ml	1.6 ml	4.8 ml	9.6 ml
Stop Solution	0.1 ml	0.8 ml	1.6 ml	4.8 ml	9.6 ml

## SUGGESTED STRIP WELL SETUP

**Table 2.** The suggested strip-well plate setup for dimethyl H4R3 quantification in a 48-assay format (in a 96-assay format, Strips 7 to 12 can be configured as Sample). The controls and samples can be measured in duplicate.

Well #	Strip 1	Strip 2	Strip 3	Strip 4	Strip 5	Strip 6
A	Blank	Blank	Sample	Sample	Sample	Sample
B	H4R3 0.5 ng	H4R3 0.5 ng	Sample	Sample	Sample	Sample
C	H4R3 1 ng	H4R3 1 ng	Sample	Sample	Sample	Sample
D	H4R3 2 ng	H4R3 2 ng	Sample	Sample	Sample	Sample
E	H4R3 5 ng	H4R3 5 ng	Sample	Sample	Sample	Sample
F	H4R3 10 ng	H4R3 10 ng	Sample	Sample	Sample	Sample
G	Sample	Sample	Sample	Sample	Sample	Sample
H	Sample	Sample	Sample	Sample	Sample	Sample

## TROUBLESHOOTING

Problem	Possible Cause	Suggestion
No signal or weak signal in both the positive control and sample wells	Reagents are added incorrectly.	Check if reagents are added in the proper order with the right amount, and if any steps in the protocol may have been omitted by mistake.
	Incubation time and temperature are incorrect.	Ensure the incubation time and temperature described in the protocol are followed correctly.
	Incorrect absorbance reading.	Check if appropriate absorbance wavelength (450 nm) is used.
	Kit was not stored or handled properly.	Ensure all components of the kit were stored at the appropriate temperature and

		the cap is tightly capped after each opening or use.
No signal or weak signal in only the standard curve wells	The standard amount is insufficiently added to the well in Step 2c.	Ensure a sufficient amount of standard is added.
	The standard is degraded due to improper storage conditions.	Follow the Shipping & Storage guidance in this User Guide for storage of <b>H4R3 Control</b> .
High background present in the blank wells	Insufficient washing of wells.	Check if washing recommendations at each step is performed according to the protocol.
	Contaminated by sample or standard.	Ensure the well is not contaminated from adding sample or standard accidentally or from using contaminated tips.
	Incubation time with <b>Diluted RDA</b> is too long.	The incubation time at Step 3d should not exceed 90 min.
	Over-development of color.	Decrease the development time in Step 4a before adding <b>SS</b> Stop Solution in Step 4b.
No signal or weak signal only in sample wells	Protein sample is not properly extracted or purified.	Ensure your protocol is suitable for histone protein extraction. For the best results, it is advised to use Epigentek's histone extraction Kit (Cat. No. OP-0006).
	Sample amount added into the wells is insufficient.	Ensure a sufficient amount of histone extracts is used as indicated in Step 2. The sample can be titrated to determine the optimal amount to use in the assay.
	Sample was not stored properly or has been stored for too long.	Ensure sample is stored in aliquots at – 80°C, with no more than 6 months histone extracts..
	Little or no di-methylH4-R3 in the sample.	This problem may be a result of many factors. If the affecting factors cannot be determined, use new or re-prepared histone extracts.
Uneven color development	Insufficient washing of the wells.	Ensure the wells are washed according to the guidance of washing and residue washing buffer is removed as much as possible.
	Delayed color development or delayed stopping of color development in the wells.	Ensure color development solution or stop solution is added sequentially and is consistent with the order you added the other reagents (e.g., from well A to well G or from well 1 to well 12).

## RELATED PRODUCTS

### Histone Extract Preparation

OP-0006      EpiQuik™ Total Histone Extraction Kit

### PRMT methyltransferase Activity/Inhibition Assay

P-3088      Epigenase™ PRMT Methyltransferase (Type II-Specific) Activity/Inhibition Assay Kit  
(Colorimetric)



110 Bi County Blvd. Ste. 122, Farmingdale, NY 11735

Tel: 1-877-374-4368 ■ Fax: 1-718-484-3956 ■ E-mail: [info@epigentek.com](mailto:info@epigentek.com) ■ Web: [www.epigentek.com](http://www.epigentek.com)

© Epigentek Group Inc. All rights reserved. Products are for research use only.